

PRESS RELEASE

Grenchen, 29th October 2018

Straight to the point: Thommen Medical presents the fully guided surgery system!

We are happy to announce the launch of the new Thommen Medical Guided Surgery Kit. Fully guided placement with Thommen Medical implants enables you to efficiently achieve your planned digital restoration, in just a few steps. Benefit from improved treatment confidence and patient comfort. The guided system will simplify the treatment process for the entire team.

Highest precision through the principle of passive guidance: Each instrument is guided precisely through the combination of integrated guidance and the proven VECTODrill™ geometry.

One drilling sleeve, perfectly matched to all cutting instruments, provides a cylindrical guide area and mechanical depth control. This integrated guidance is significantly easier in comparison to indirectly guided systems (e.g. tray systems). Additionally, the same surgical cassette is utilized for everything: surgical treatment, processing/sterilization and storage.

Contact

Thommen Medical AG
Neckarsulmstrasse 28
CH-2540 Grenchen
Tel. +41 61 965 90 20
Fax +41 61 965 90 21
presse@thommenmedical.com
www.thommenmedical.com

The press release is also available online:

<http://www.thommenmedical.com/en/company/media/press-releases.html>

About Thommen Medical

Thommen Medical is a Swiss designer and manufacturer of a dental implant system known for incredible precision and impactful design. We are unique with an incredible heritage of over 30 years of Swiss precision, innovation, and functional design that has formed the basis for our products. Of high quality yet simple to use, our products deliver predictable, aesthetic results without compromise. Thommen Medical's headquarters and private manufacturing facility are in Grenchen, Switzerland, a location that is the heart of Swiss precision manufacturing. Thommen Medical is an international company with subsidiaries and distributors throughout the world.