

GUIDED SURGERY with steco

The StecoGuide double sleeve system for Thommen Medical

- StecoGuide double sleeve system for VECTOdrill™ twist drills stainless steel
 - Guided drill sequence up to VECTOdrill™ twist drill Ø 3.5 mm
 - Flexible placement of the sleeve (offset)
 - Depth control via lasermarking
 - Certified for the following planning systems: 3Shape “Implant Studio”, Swissmeda “SMOP”
-

The StecoGuide sleeves for the Thommen Medical VECTOdrill™ instruments are directly available via steco:

steco-system-technik GmbH & Co. KG
Kollastr. 6
22529 Hamburg
Germany
Tel: +49 (0)40 55 77 81-0
Fax: +49 (0)40 55 77 81-99
E-Mail: info@steco.de
www.steco.de




And other local distributors




SYSTEM OVERVIEW

Double sleeve system for drilling sequences



Pilot sleeve for pilot drilling only



- 1 Inner sleeve M.27.25.D200L6
- 2 Inner sleeve M.27.25.D280L6
- 3 Outer sleeve M.27.25.D350L6



- I VECTOdrill™ pilot drill Ø 2.0 mm
- II VECTOdrill™ twist drill Ø 2.8 mm
- III VECTOdrill™ twist drill Ø 3.5 mm

- 1 Inner sleeve M.27.25 D200L6

- I VECTOdrill™ pilot drill Ø 2.0 mm



SWISS CRAFTSMANSHIP FOR DENTAL PROFESSIONALS

Thommen Medical AG | Headquarters

Neckarsulmstrasse 28 | CH-2540 Grenchen | Tel. +41 61 965 90 20 | Fax +41 61 965 90 21 | info@thommenmedical.com
www.thommenmedical.com